

AquaStar

Caloric Irrigator



- Ø QUIET
- Ø PRECISE
- Ø COMPACT
- Ø SAFE
- Ø DEPENDABLE

The AquaStar water caloric irrigator delivers a precisely controlled temperature and flow rate of cool or warm water for vestibular stimulation during caloric testing. In less than a minute from turning on the AquaStar the temperature stabilizes at the desired setting as indicated on the digital display. A timing mechanism then provides audible and visual cues for irrigating.

The AquaStar connects to the water supply in your exam room eliminating noisy pumps and the need to refill external tanks. Preset temperatures of 30°C and 44°C are switch selectable. AquaStar is designed with built-in filters to ensure a clean water stream, replaceable silicone delivery tips and with thermal and electrical protective circuitry for ultimate patient safety.

Specifications

Temperature precision	± 0.5°C	Temperature Increase Stabilize Time	45 seconds
Temperature Settings	1°C above inlet temp - 44°C	Temperature Decrease Stabilize Time	45 seconds
Readout	Digital	Case Dimensions	12w x 4.5h x 12d (inches)
Flow rate (adjustable)	350 mL / minute	Weight	12 pounds
Timer Range	20 to 60 seconds	Power Requirements	110VAC, 12A
		Also Available in	220VAC, 8A



10 Kemp Drive
Chatham, IL 62629
800.334.4154 U.S. and Canada
217.483.4533 FAX
www.micromedical.com

©2002, Micromedical Technologies, Inc.
02/02 r08/02 r08/06 Printed in U.S.A.